


<b>Date</b>	06 Jun 2012	<b>Service Note #</b>	4300-1048	<b>Revision</b>	06 Jun 2012
<b>Product</b>	Model 4300	<b>Description</b>	2500_7100 Bath System Configuration		
<b>Created By</b>	J. Yangco				

<b>Release</b>	<input checked="" type="checkbox"/>	<b>Internal</b>	<input checked="" type="checkbox"/>	<b>Distributors</b>	<input checked="" type="checkbox"/>	<b>Customers</b>
----------------	-------------------------------------	-----------------	-------------------------------------	---------------------	-------------------------------------	------------------

<b>Parts Required</b>	<b>Serial Numbers Affected</b>
4300 Firmware Revision 3.02 or later	N/A

1. Connect the dissolution bath RS-485 I/O port (using 2700-0162A cable) in to the 4300 pump RS-485 port as shown below.
2. Change the instrument(s) address (refer to the instrument manual). If you have two dissolution baths change the address to 03 for Bath A and 04 for Bath B.
3. Log in to the 4300 as part of the Manager group.
4. On the main menu, select "ADMIN" then "CONFIG". The "CONFIGURATION SETTINGS" screen will be displayed.
5. Scroll thru the screen so you can change the address for BATH A to 03. Do the same if you have a second bath (BATH B set to 04).
6. Change the bath TYPE to reflect the dissolution bath model (2500 or 7100) that is connected to the 4300.



**Make sure that the firmware on the 4300 Autosampler has Revision 3.02.**