



6	5510-0144	SS DOWEL PIN 1/8" X 3/8"	1
5	5510-0146	SS CUP PT SET SCREW #4-40 X 1/4"	1
4	5510-0105	S/S SOC HD CAP SCREW #6-32 X 3/8" LG	4
3	5510-0133	NYLON THUMB SCREW 6-32 X 3/8"	1
2	3200-0241	DETECTOR INLET BLOCK	1
1	3200-0240	DETECTOR INLET BASE	1
ITEM	PART NO.	DESCRIPTION	QT'Y

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO OTHERS FOR ANY PURPOSE NOR USED FOR MANUFACTURING PURPOSES WITHOUT WRITTEN PERMISSION FROM DISTEK.

DISTEK

121 NO. CENTER DR.
NORTH BRUNSWICK, N.J. 08902

TITLE

DETECTOR INLET ASSEMBLY

OPT-DISS 410

SCALE

3/4

SIZE
A

SHEET OF
11

REV
A

DRAWING NO.

2800-0496