



6	6500-0012	1 LITER WIDE MOUTH BOTTLE, POLYCARBONATE	1
5	5720-0260	PVC SUCTION LINE STRAINER, 40 MESH, 3/8 NPT	1
4	3250-0165	FILTER EXTENSION PIPE	1
3	3250-0120	GASKET, FILTER BOTTLE	1
2	3250-0164	FILTER HEAD	1
1	5720-9021	ELBOW FITTING, 3/8 NPT, 5/8 HOSE ID BARB	2
ITEM	PART NO.	DESCRIPTION	QTY

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO OTHERS FOR ANY PURPOSE NOR USED FOR MANUFACTURING PURPOSES WITHOUT WRITTEN PERMISSION FROM DISTEK.

**DISTEK**

121 NO. CENTER DR., NORTH BRUNSWICK, N.J. 08902

TITLE

FILTER\_ASSEMBLY  
VIP\_4400

SCALE

NTS

SIZE

A

SHEET

1

OF

1

REV.

A

DRAWING NO.

REVIEW\_ONLY