



4851-0007
LCD DISPLAY

5	3200-0229	SPRING CAP, DISPLAY ARM	1
4	5590-0015	S/S SPRG .112W.D. .609OD .63LG	1
3	5510-0108	S/S PFH SCW #10-32 X 1 1/4" LG	1
2	2800-0204	LCD ASSEMBLY	1
1	3300-0003_M	DISPLAY ARM	1
ITEM	PART NO.	DESCRIPTION	QTY

CONFIDENTIAL INFORMATION
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION
AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO
OTHERS FOR ANY PURPOSE NOR USED FOR
MANUFACTURING PURPOSES WITHOUT WRITTEN
PERMISSION FROM DISTEK.

DISTEK
121 NO. CENTER DR.
NORTH BRUNSWICK, N.J. 08902

TITLE		SCALE	
DISPLAY ARM ASSEMBLY SYMPHONY 7100/7300		5/12	
		SIZE B	SHEET 1 OF 1 REV. A
		REVIEW ONLY	