


A B C D E F G H J

1  
2  
3  
4  
5



ITEM	PART NO.	DESCRIPTION	QTY
7	3300-0001-M	VESSEL PLATE, MACHINED	1
6	3200-0220	TOP PLATE	1
5	5300-0007	TIMING BELT, .080P, 8.0L, .25WIDE	1
4	2800-0193	MOTOR-PULLEY ASSEMBLY	1
3	3200-0213	GUIDE ROD	2
2	3250-0216	CABLE GUIDE	1
1	2800-0195	BRAKE ASSEMBLY	1

CONFIDENTIAL INFORMATION  
 THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION  
 AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO  
 OTHERS FOR ANY PURPOSE NOR USED FOR  
 MANUFACTURING PURPOSES WITHOUT WRITTEN  
 PERMISSION FROM DISTEK.

**DISTEK**  
 121 NO. CENTER DR.  
 NORTH BRUNSWICK, N. J. 08902

TITLE  DRIVE SYSTEM ASSEMBLY  SYMPHONY 7100/7300	SCALE 1/4		
	SIZE B	SHEET 1 OF 1	REV C
	REVIEW ONLY		