



3	3200-0217	MOTOR BRACKET	1
2	5400-0001	PLASTIC PULLEY, 15T	1
1	5100-0015	BRUSHLESS DC MOTOR	1
ITEM	PART NO.	DESCRIPTION	QT'Y

CONFIDENTIAL INFORMATION  
 THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION  
 AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO  
 OTHERS FOR ANY PURPOSE NOR USED FOR  
 MANUFACTURING PURPOSES WITHOUT WRITTEN  
 PERMISSION FROM DISTEK.

**DISTEK**

121 NO. CENTER DR.  
 NORTH BRUNSWICK, N.J. 08902

TITLE  
 MOTOR-PULLEY ASSEMBLY  
 SYMPHONY 7100/7300

SCALE 1/2  
 SIZE A SHEET 1 OF 1 REV C

REVIEW ONLY