



ITEM	PART NO.	DESCRIPTION	QTY
4	5535-0006	#4-40 KEPS NUT, S.S	4
3	5510-1145	FLANGED BUTT. HD SCREW, #4-40 x 3/8" LG	4
2	5500-0002	BALL STUD RECEIVER, 12 LB. PULL OUT	2
1	3350-0012	FRONT COVER BOTTOM	1

CONFIDENTIAL INFORMATION
 THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION
 AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO
 OTHERS FOR ANY PURPOSE NOR USED FOR
 MANUFACTURING PURPOSES WITHOUT WRITTEN
 PERMISSION FROM DISTEK.

DISTEK
 121 NO. CENTER DR.
 NORTH BRUNSWICK, N.J. 08902

TITLE		SCALE	
FRONT COVER BOTTOM ASS'Y ECLIPSE 5300		1/2	
		SIZE B	SHEET OF 11
		REV	A
		REVIEW ONLY	