


A B C D E F G H J

1  
2  
3  
4  
5



USE LOCTITE 425 ON ITEM 2 & 3 THREADS

LEFT SIDE COVER ASSEMBLY 2800-0423		
ITEM	PART NUMBER	REPLACE WITH
1	3250-0297	3250-0298

3	5500-0003	BALL STUD, FLAT HEAD, MALE THREAD	1
2	5500-0001	BALL STUD, MALE THREAD	5
1	3250-0297	SIDE COVER RIGHT	1
ITEM	PART NO.	DESCRIPTION	QTY

CONFIDENTIAL INFORMATION  
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION  
AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO  
OTHERS FOR ANY PURPOSE NOR USED FOR  
MANUFACTURING PURPOSES WITHOUT WRITTEN  
PERMISSION FROM DISTEK.

**DISTEK**

121 NO. CENTER DR.  
NORTH BRUNSWICK, N.J. 08902

TITLE  
**SIDE COVER, RIGHT (NOT SHOWN)**  
**ECLIPSE 5300**

SCALE **1/3**  
SIZE **B** SHEET OF **11** REV **B**  
**REVIEW ONLY**