



ITEM	PART NO.	DESCRIPTION	QTY
9	3200-0306	FRAME PLATE	2
8	3100-0216	BRACKET, SLOTTED, STAGE HORZ SENSOR	1
7	3100-0215	INTERNAL COVER	1
6	3100-0034	BRACKET, END, STAGE HORZ. SENSOR	1
5	2700-0054	SENSOR, HALL, FRONT COVER INTERLOCK #3	1
4	2700-0029	WIRING HARNESS, TUBE SENSOR, NEEDLE LED	1
3	2700-0028	WIRING HARNESS, TUBE LED, NEEDLE SENSOR	1
2	2700-0026	SENSOR, RACK PRESENCE	1
1	2700-0009	SENSOR, HALL, SIDE COVER INTERLOCK #2	2

CONFIDENTIAL INFORMATION  
 THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION  
 AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO  
 OTHERS FOR ANY PURPOSE NOR USED FOR  
 MANUFACTURING PURPOSES WITHOUT WRITTEN  
 PERMISSION FROM DISTEK.

**DISTEK**  
 121 NO. CENTER DR.  
 NORTH BRUNSWICK, N.J. 08902

TITLE		SCALE	
SUPPORT STRUCTURE ASS'Y ECLIPSE 5300		1/2	
		SIZE B	SHEET OF 11 REV B
		REVIEW ONLY	