



ITEM	PART NO.	DESCRIPTION	QTY
6	5720-9020	COMPRESSION FITTING	1
5	5720-7076	HEATER SENSOR NUT	1
4	5720-0047	INSTANT FITTING ELBOW	1
3	4921-1008	HEATER SENSOR WIRE ASSY	1
2	3200-0270	HEATER BLOCK	1
1	2700-0313	HEATER CART ASM 110V	1

CONFIDENTIAL INFORMATION
 THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION
 AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO
 OTHERS FOR ANY PURPOSE NOR USED FOR
 MANUFACTURING PURPOSES WITHOUT WRITTEN
 PERMISSION FROM DISTEK.

DISTEK
 121 NO. CENTER DR.
 NORTH BRUNSWICK, N.J. 08902

TITLE		SCALE	
HEATER BLOCK ASSY 110V TCS 0500		2/3	
		SIZE B	SHEET OF 1 REV. A
		REVIEW ONLY	