

A B C D E F G H J


1

2

3

4

5



2851-0024, 230 VAC UNIT		
ITEM	PART NUMBER	REPLACE WITH
1	3500-0345	3500-0347

FINISHED ASSEMBLY PART NUMBER: 2851-0015

7	5510-1084	S/S FHP SCREW #6-32, 7/16 LG	4
6	3350-0003	HEATER RIVET	8
5	3251-1080	SPRING COVER	4"
4	3251-1109	EXTENSION SPRING	2
3	5535-2108	S/S HEX NUT #6-32 WITH NYLON INSERT	4
2	3351-1064	ENCLOSURE, FLUTED HEATER	1
1	3500-0345	FLEXIBLE HEATER, 110V	1
ITEM	PART NO.	DESCRIPTION	QTY

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION AND SUCH INFORMATION IS NOT TO BE DISCLOSED TO OTHERS FOR ANY PURPOSE NOR USED FOR MANUFACTURING PURPOSES WITHOUT WRITTEN PERMISSION FROM DISTEK.

DISTEK

121 NO. CENTER DR., NORTH BRUNSWICK, N.J. 08902

TITLE	SCALE	NTS
HEATER_ASSY_110V 6100_7100_2500	SIZE B	SHEET 1 OF 1 REV E
	DRAWING NO.	REVIEW_ONLY