



Release Note

Date	19 Jan 2011		Release Note #	6x00-1001	ECO Ref #	0861
Product	Evolution 6x00		Description	Stage Board Firmware 2.15		
Revision Change	From	2.14	To	2.15	Created By	JSB

Release	<input checked="" type="checkbox"/>	Internal	<input checked="" type="checkbox"/>	Distributors	<input type="checkbox"/>	Customers
----------------	-------------------------------------	-----------------	-------------------------------------	---------------------	--------------------------	------------------

Assembly / Part #	N/A (Firmware Revision)
Serial #	Evolution 6100 (Bathless) – 6102011 or Later Evolution 6300 (Water Bath) – 6301374 or Later
Description	<p>Due to a change in the LCD controller the Stage Board Firmware was updated to be compatible with the new LCD. Certain graphics, text and initialization codes were modified.</p> <p>This Firmware is only compatible with the new version of the LCD (Distek P/N 4851-0006) and is not backwards compatible with the obsolete LCD (Distek P/N 4851-0003).</p>