

Date	06 Sept 2011		Release Note #	4300-1012	ECO Ref #	0883
Product	Evolution 4300		Description	EVO 4300 FW 3.01		
Revision Change	From	3.00	To	3.01	Created By	J. Brinker

Release	<input checked="" type="checkbox"/>	Internal	<input checked="" type="checkbox"/>	Distributors	<input checked="" type="checkbox"/>	Customers
----------------	-------------------------------------	-----------------	-------------------------------------	---------------------	-------------------------------------	------------------

Assembly / Part #	N/A
Serial #	4301815 or Later
Description	Provisions were made to the Evolution 4300 firmware to allow control of Distek's Symphony Dissolution Systems. Instrument types were added to the corresponding pull down menus in the user interface.



Currently the symphony Dissolution System in conjunction with the Evolution 4300 is only supported for a maximum of 7 dissolution vessel positions. Support for 8 dissolution vessel positions will be addressed in a future Evolution 4300 Firmware release.