









Date	23 Apr 2014	Service Note #	6x00-2002	Updated	NA
Product	Evolution 6100 / 6300		Created By	J. Brinker	
Description	6x00 LCD Display and Stageboard PCB FW Compatibility Matrix				

Release	<input checked="" type="checkbox"/>	Internal	<input checked="" type="checkbox"/>	Distributors	<input checked="" type="checkbox"/>	Customers
----------------	-------------------------------------	-----------------	-------------------------------------	---------------------	-------------------------------------	------------------

Parts Required	Serial Numbers Affected
N/A	All

RTD PCB Firmware, Stageboard EPROM and Stageboard PLD Firmware - Compatibility Matrix

The RTD circuit board (2451-5108) was modified due to a component change which required a Firmware revision on the Stageboard PLC chip and EPROM for compatibility. The table below lists the compatibility matrix.

	RTD PCB Rev 1.0	RTD PCB Rev 1.1 or Higher
Stageboard PLD Firmware Rev 1.0		
Stageboard PLD Firmware Rev 1.1		
Stageboard EPROM Firmware Rev <= 2.13		
Stageboard EPROM Firmware = Rev 2.14 or Higher		

Stageboard 2400-5154



RTD Board 2451-5108

